

Gebrauchsinformation

Happy Taq



IVD

REF 70976

- Inhalt:** 50 µl Taq Polymerase, gebrauchsfertig
250 units (5 units/µl)
- Spezifikation:** Rekombinantes Enzym
Thermostabile 5' → 3' DNA-Polymerase
5' → 3' Exonuclease Aktivität
Keine 3' → 5' Exonuclease Aktivität
Enzym fügt Nukleotide (fast ausschließlich Adenosin) an das 3'-Ende der DNA
- Lagerung:** ≤ -20°C
- Verwendung:** Für den Einsatz mit HISTO TYPE Kits der BAG Health Care GmbH.
Für die Testdurchführung und weitere Informationen siehe Gebrauchsinformation für die HISTO TYPE SSP Kits.

Erklärung der Symbole auf den Etiketten			
	Lagertemperatur	IVD	In-vitro-Diagnostikum
	Verwendbar bis	REF	Bestell-Nr.
LOT	Lot-Nr.	RTU	Gebrauchsfertig
		Gebrauchsinformation beachten	
TAQ POLYMERASE		Taq-Polymerase	

Version: 1/2014 / Stand: 2014-05

Instructions for use

Happy Taq



IVD

REF 70976

- Content:** 50 µl Taq Polymerase, ready to use
250 units (5 units/µl)
- Specification:** Recombinant enzyme
Thermostable 5' → 3' DNA polymerase
5' → 3' exonuclease activity
No 3' → 5' exonuclease activity
Enzyme adds nucleotides at the 3'-end of the DNA (almost exclusively adenosine)
- Storage:** ≤ -20°C
- Application:** For use with HISTO TYPE kits of BAG Health Care GmbH.
For test procedure and further information see instructions for use for HISTO TYPE SSP kits.

Explanation of symbols used on labelling			
	Storage temperature	IVD	For in vitro diagnostic use
	Use by	REF	Catalogue number
LOT	Batch code	RTU	Ready to use
		Consult instructions for use	
TAQ POLYMERASE		Taq-Polymerase	

Version: 1/2014 / Issue: 2014-05

Istruzioni d'uso

Happy Taq



IVD

REF 70976

- Contenuto:** 50 µl Taq Polimerasi, pronta all'uso
250 unità (5 unità/µl)
- Specifiche:** Enzima ricombinante
Termostabile 5' → 3' DNA polimerasi
Attività 5' → 3' esonucleasica
No attività 3' → 5' esonucleasica
L'enzima aggiunge nucleotidi all' estremità 3' del DNA (quasi esclusivamente adenosina)
- Conservazione:** ≤ -20°C
- Applicazioni:** Per l'uso con i kit HISTO TYPE della BAG Health Care GmbH.
Per la procedura del test e ulteriori informazioni, vedere istruzioni d'uso per i kit di HISTO TYPE SSP.

Spiegazione dei simboli presenti sulle etichette			
	Temperatura di conservazione	IVD	Per uso diagnostico in vitro
	Utilizzare fino a	REF	Codice prodotto
LOT	Codice lotto	RTU	Pronto per l'uso
		Consultare le istruzioni d'uso	
TAQ POLYMERASE		Taq-Polimerasi	

Versione: 1/2014 / Data di rilascio: 2014-05

BAG Health Care GmbH

Amtsgerichtsstraße 1-5 Tel.: +49 (0) 6404 / 925-0 www.bag-healthcare.com
35423 Lich / Germany Fax: +49 (0) 6404 / 925-250 info@bag-healthcare.com

Auftragsannahme/Ordering:

Tel.: +49 (0) 6404 / 925-450
Fax: +49 (0) 6404 / 925-460
verkauf@bag-healthcare.com

Customer Service:

Tel.: +49 (0) 6404 / 925-125
Fax: +49 (0) 6404 / 925-421
service@bag-healthcare.com